

BIOLOGICAL SAFETY CABINET

Protect Laboratory Workers and Surrounding Environment



BIOBASE BSC-1100IIB2-BIOSAFETY CABINET 700 CLASS 11 A2

- External Size: (W*D*H): 700*650*1230mm (without base stand)
- Max Opening: 370mm (15")
- HEPA Filter: Two, 99.9995% efficiency at 0.3µm filter life indicator.
- Front window: Motorized. Single layer toughened glass 5mm Anti UV
- Display: LCD display
- Visual and audio alarm: abnormal airflow velocity, filter replacement, window over height, high filter pressure alarm.

BIOBASE BSC-1800IIB2-PRO BIOSAFETY CABINET 1800 CLASS II B2

- External Size (W*D*H): 1873*780*2310mm
- Filter life indicator
- Front window: Motorized. Two layer laminated toughened glass.
- Display: 7 inch touch screen display: exhaust filter and downflow filter pressure filter and UV lamp working time, inflow and downflow velocity, filter life, humidity and temperature, system working time etc.
- UV lamp: UV timer, UV life indicators, emissions of 253.7 nanometers for most efficient decontamination.
- Visual and audio alarm: abnormal airflow velocity, filter replacement, window over height, high filter pressure alarm.



BIOBASE BSC-1100IIB2-BIOSAFETY CABINET 700 CLASS 11 A2

- External Size: (W*D*H): 700*650*1230mm (without base stand)
- Max Opening: 370mm (15")
- HEPA Filter: Two, 99.9995% efficiency at 0.3µm filter life indicator.
- Front window: Motorized. Single layer toughened glass 5mm Anti UV
- Display: LCD display
- Visual and audio alarm: abnormal airflow velocity, filter replacement, window over height, high filter pressure alarm.

To Request for Quote:

WhatsApp: +254 739 643897