

Cold Chain Equipment



Apex Vaccine Carrier Box AIVC-44

- Apex Vaccine Carrier Box AIVC 44



BLOOD BANK 398L FRIDGE

- Voltage/Frequency 220/50Hz
- Specification: 1330*740*2120mm
- Ambient Temperature: 10~32°C
- Digital Display



Classic Cold Pack & Buffer Board

- Cold pack for the CoolMeds Classic 2 – 8°C Medication Case
- Each cold pack comes with a free buffer board
- Made from long-lasting coolant
- Keeps medication cool at 2-8°C when used with the CoolMeds Classic case
- Stays cold for 12-18 hrs
- x 1 cold pack and x 1 buffer board

To Request for Quote:

WhatsApp: +254 739 643897



COLD BOX

- Range: Long
- Vaccine Storage Capacity: 22.5 Litres
- Weight Fully Loaded: 48.2 Kgs
- Weight Empty (with empty ice pack): 23.55 Kgs
- External dimensions(WxDxH)*: 76.5 X 61.2 X 5.5cms
- Internal dimensions(WxDxH)*: 51.8 X 37.0 X 27.7cms
- Vaccine Storage dimensions (WxDxH): 45.5 X 31.0 X 16.0cms



COOLER/COLD BOX 14L

- Suitable for: Blood center, testing center, vaccine center, bio
- Temp Range: 2~8°C
- Cooling Time: 3-5h



COOL PACK 255×150 mm

Available in two size: Regular and Large

Applications:

- Control high-temperature fever
- Instant relief in sports injuries
- A soothing effect on bruises and bumps
- A relaxant for eye strain and headaches
- Effective on sprains
- Quick relief from toothaches
- Also recommended in hangover, insect bites, stings, sun stroke tension and migraine



DIGITAL THERMOMETER AND HYGROMETERS

- Temperature measurement: 58°F-158°F/-50°C-70°C;
- Humidity Measurement: 10% to 99% RH; Accuracy of temperature: 0.1°C; accuracy of humidity:1%RH
- 5 Function: VADIV temperature thermometer displays humidity percentage, clock, alarm, indoor and outdoor temperature, powered by a 1.5v AAA battery (NOT included)
- Convenient Design: You can hang our electronic temperature humidity gauge or place it on table by its standing and wall hanging design
- Dimension: (approx) W10 x H10x D2 cm
- Waterproof Probe: Comes with 1M sensor cable for outside temperature, so the weather thermometer hygrometer is ideal for aquarium, school, greenhouse, garage, warehouse, baby room, laboratory, garden



FRIDGE TEMP RECORDERS WITH USB

- Power LCD indication
- Software-free Online configuration
- Multiple Sensors for Temperature, humidity and illumination



FRIDGE THERMOMETERS

- Overall dimension: 91mm (length)
- 20.8mm (width) 10.6mm (thickness)
- Temperature range: -30-40
- Measurement accuracy: +/-2oc
- Method of use: hanging or adsorption



FRIDGE THERMOMETERS (HANGING OR ADSOPTION)

- Temperature range: -30~30°C – 20~80°F
- Accuracy: +/- 2°
- Material: Stainless steel case
- Display: °C/°F display



PHARMACY FRIDGE 138L

- **Volume:** 138L
- **Inner Temp:** 2~8°C
- **Voltage:** 220V/50Hz
- **Alarm type:** Sound/Light