



Dental Equipment



AMALGAM CARRIER

- Dental Amalgam Filling Carrier Regular/Jumbo 2mm/3mm



APEX LOCATOR

-Apex Locator WOODPEX III PRO



BONE CHISELS

To Request for Quote:

WhatsApp: +254 739 643897



BURNISHER

- Dental Ball Burnisher



CORDLESS GUTTA PERCHA OBTURATION

- Cordless Gutta Percha Obturation
- C FILL Obturation Gun



COXO PROPHY CONTRA ANGLE

- Push button, Max 40,000rpm, CA bur 2.35mm




DENTAL AIR COMPRESSOR 60L

- "Power: 600W (0.75hp)
- Flow: 118L/min at 0 bar (4.2cfm)
- Volt./Hz: 110-240V, 50-60Hz
- Noise Level: 52dB
- Pressure: 5-8bar
- Tank: 60L (15.85 gallon)
- G.W.: 32kg, N.W.: 24kg
- Dimension: 410*410*550mm"



DENTAL EXPLORERS

- Made of High quality Medical grade Stainless Steel.
- Ergonomic design, double-ended tips, single-ended tips with hollow handle helps reduce hand fatigue



DENTAL LABORATORY POLISHING LATHE

- Dental Laboratory Polishing Lathe



DENTAL TWEEZERS

- Single fold with pin
- Made of 304 stainless steel, which is very strong and rustproof and excellent flexibility



DENTAL UNIT COMPLETE

- "Dental LED light
- Ceramic spittoon
- Auto spittoon flush and cup filler sensor control system
- Assistant control box (13-15) with 2 tubes
- 3 way syringe (cold/hot water) 2 pcs
- Saliva ejector with high and low suction
- 24V noiseless DC motor chair with one arm
- Instrument tray with air brake with 4 tubes
- Doctor chair
- Round foot pedal (04-01)
- X-film viewer



HAEMOSTATIC FORCEPS

- Haemostatic Forceps



IMPRESSION TRAY ALUMINIUM

- Impression Tray Aluminium



LUXATING ELEVATORS

- Luxating Elevators



MATRIX BAND

- Dental Matrix Bands Strips



MOUTH MIRRORS

- Mouth Mirrors for precise dental examinations, High-quality stainless steel construction ensures clarity and durability.



NEEDLE HOLDERS

- Needle Holders



PLASTIC FILLING INSTRUMENTS

- Plastic Filling instruments



WAX CARVERS

- Wax Carvers