

General Ward Equipment



ABS MEDICINE TROLLEY

- **Specification:** 750*475*930mm
- **Material:** Stainless Steel
- **Wheel Type:** Swivel Casters
- **Type:** Treatment Cart
- **Layer:** 5 Layers
- **Feature:** Easy cleaning



BED PAN

- Perfection Type – with cover, ss



DEHUMIDIFIER 138L

- Microcomputer control, intelligent operation
- LED display, indicate humidity
- Timing range: 1~24 hours
- Pipe continue drainage function
- Low temperature available, efficient automatic defrosting system

To Request for Quote:

WhatsApp: +254 739 643897



DRIP STAND

- Robust construction: sturdy frame providing stability for secure support of IV fluid bags.
- Height adjustability: easily adaptable to meet varying patient requirements and clinical settings
- Nobility with lockable wheels: convenient movement and stability ensured by wheels with locking capability.
- Drip rate control (select models): allows precise adjustment for accurate medication administration.
- Hygienic and standardized: built with corrosion-resistant materials, it is a standard tool in medical settings for intravenous treatments.



ELECTRIC NURSING BED

- Four splits side rail, drop down under the bed platform
- Hand remote control
- Detachable head/footboard with easy lock
- With articulated ABS soft joint between different platform sections.
- Central brake system with a wheel diameter of 125mm.
- Auto-counter to provide the patient from sliding onto foot end



EXAMINATION TABLE, ELECTRIC (TWO SECTIONS)

- Approx. Size : 190L x 64W x 48-81H cm
- Height and backrest Adjustable through electric remote control.
- Two section foam padded, upholstered top
- Pre treated & Epoxy powder coated
- Mild Steel Framework



GYNECOLOGICAL TABLE, ELECTRIC

- Length 140cm
- Width 60cm
- Minimum Height 60cm
- Maximum Height 100cm
- Back Lift : Electric Control
- Up/Down-Height Adjustment : Electric Control
- Framework : Mild Steel
- Finish : Beautiful Epoxy Powder Coating
- Castors : 150 mm dia
- Brakes : Individual braking Castors
- Control Panel : Foot operated controls for easy adjustments
- Hand grip : Two hand grips provide good support for the patient



HTC-01 TOILET CHAIR

- Easy to Clean & Compact Storage with Lock for added security.
- This Commode Easily opens & folds and will fold flat for convenient storage and transportation.



INFUSION PUMP

- Double CPU, safer infusion
- Brightness LCD displays visually
- Numeric keypad, user friendly
- 10 level occlusion pressure detection adjustable
- Occlusion air-in-line, audio visual alarm volume: 10 settings
- 3000Ah Lithium battery, backup life: ≥10 hours
- KVO adjustable 0.1-5ml/h



JAY-10 OXYGEN CONCENTRATOR SINGLE FLOW

- Flow Rate: 10 LPM
- Oxygen Concentration: 95%
- Sound Level: <52db
- Dimensions: 45x65x75



LABOUR & DELIVERY ROOM(LDR) BED, ELECTRIC

- Size : 188L x 90W x 60-80H cm (when fully extended)
- Size:150L x 90W x 60-80H cm (when leg section under the middle section)
- Functions : Five Functions
- Back Lift : Electric Control, 0 to 75 degree
- Trendelenburg : Electric Control, 0 to 12 degree
- Up/Down-Height Adjustment : Electric Control
- Head & Foot Panels : Elegant High Quality Moulded Panels.
- Side Support : Moulded safety side panels.
- Top : Three Sections
- Leg section can be completely recessed under the body section.
- Framework : Mild Steel.
- Castors : 150 mm dia with central locking mechanism



MANUAL DIALYSIS BED

- All-steel frame provides strong bearing capability, safe working load up to 240kg
- High density polyurethane foam for mattress with moderate elasticity and durability provides additional comfort for patients having to sit for long periods.
- Soft upholstery made of breathable PU material, corrosion, easily clean and disinfection.
- Three upholstery colors to choose from, warm and light colors of upholstery contributes to patients' mental well-being.
- Mute medical castors of diameter 12.5 cm with central brakes.
- The detachable headrest pillow can be adjusted according to the height of each patient.
- The collapsible aluminum guardrails and armrests are easy to operate with one button



MANUAL SINGLE CRANK BED WITHOUT MATTRESS

- Back angle: 0°-80°
- Strip steel bed frame (powder coating surface treatment)
- Detachable PP blow molded head & foot board with anti-bump corners.
- Aluminum guardrail of 5 bars.
- 5" diameter noiseless casters with independent lock
- ABS crank handle
- Safety load: 200kg



URINAL / URINE POT – FEMALE, SS



WALL MONITOR HOLDER WITH BASKET

- High quality
- Long lasting
- Made in Aluminum extrusion
- Movement: Up and down, left and right turn
- Meet most of the brand
- Monitor stand (adjustable with rolling base)
- Integrated cable management
- One quick release baskets
- Single grip handle
- Seamless design small-medium/medium-large