



Laboratory Equipment



BLOOD TUBE ROLLER ORDINARY

- Speed: 8-32rpm
- Capacity: 1-15pcs/time
- Diameter of tube: $\varnothing 13 \pm 0.1 \text{mm}$
- Diameter of plate: $\varnothing 220 \text{mm}$



BLOODTUBE ROLLER-360

- Operation mode: 360 reverse rotation
- Tilt angle of the plate : 38°
- Speed range: 6-32rpm
- Clip Qty: Large 12pcs/Medium 24pcs/small:24pcs
- Running mode: Timing/continuous
- Plate diameter: 300mm



CENTRIFUGE 6T ROUND 6*20ml

- Power Supply: 110V 60HZ
- Maximum Speed: 4000r/min
- Capacity: 20ml*6 (total six tubes)
- Timer Range: 0-60min

To Request for Quote:

WhatsApp: +254 739 643897



CLOT ACTIVATOR TUBE

- Clot Activator Tube. 13*75mm, Glass



DYMIND DH 36 AUTOMATED HEMATOLOGY ANALYZER

- Product DYMIND Hematology Analyzer 3
- Model: DYMIND Dh36
- After-sale service: Online technical
- Storage: up to 50000 records
- Function: Blood Analysis System
- Number of parameters: 21



ESR ANALYZER

- Large color display, touch screen
- 10 channels, up to 80T/H
- Two working method: 30 or 60 mins
- ESR curve display and print out
- Barcode reader optional



GETEIN 1160 IMMUNOFLUORESCENCE QUANTITATIVE ANALYZER

- Cardiac Markers: cTnI, TnT, BNP, NT-proBNP, H-FABP, CK-MB/cTnI, CK-MB/cTnI/Myo, CK-MB/cTnI/H-FABP Coagulation Markers: D-Dimer
- Inflammation: hs-CRP+CRP, PCT, SAA, IL-6
- Renal function: CysC, mAlb, NGAL, B2-MG
- Diabetes: HbA1c
- Metabolic Markers: tPSA, fPSA, AFP, CEA
- Reproduction/Fertility: HCG+β, FSH, LH, PRL, AMH, Prog, Testosterone, E2
- Infectious Disease: Anti-HCV, Anti-TP, Anti-HIV, SARS-Cov-2 Neutralizing Antibody, HBsAg, Anti-HBs, SARS-Cov-2 Antigen, H-Pylori, Influenza A/B, Dengue NS1 Ag



LABORATORY INCUBATOR 30L

- Power supply: 220V 50HZ 110V 60HZ 1
- Temperature range: RT+5~65°C
- Inner dimension: 300*300*350mm



LABORATORY STOOL WITH BACKREST

- Anti-static,
- Adjustable height,
- Stainless steel



LABORATORY TIMERS

- Material: Plastic
- Type: Digital



MICROPIPETTES HOLDER/STAND

- Material: Plastic



MINDRAY BC 5000 AUTOMATED HEMATOLOGY ANALYZER

- Operation Mode: fully automatic
- Differential Type: 5-part
- Model Name/Number: Bc5000
- Display Type: 10.4 inch TFT Touch screen
- Throughput: 40 samples per hour
- Principles: Impedance method for RBC and PLT counting
- Reagent Diluent, DIFF lyse, LH Lyse, Probe Cleanser
- Dimension depth (400mm x Width (320mm) x Height (410mm)
- Weight: 24kg
- Voltage: 100v – 240V
- Frequency: 50Hz/60Hz
- Air Pressure: 70kPa~106kPa
- humidity: 20%~85%



PAPSMEAR KIT

- The pap smear kit is used for efficient and adequate collection of uterine cervical specimens
- The enclosure package promotes specimen safety during collection and transportation
- Contents: Cyto brush/ Endocervical brush, Cotton applicator/ Cotton stick/ Cotton tipped applicator, Wooden cervical scraper, Glass slide
- Clear & sure Pap Smear Kit contains one wooden spatula, cervical brush & Swab stick sterile in one blister pack.
- Pap smear Test Kit uses for women



URISTIK (100 STRIPS)

- PH, Nitrite, Protein, Glucose, Occult Blood, Bilirubin, Urobilinogen, Ketone, Specific Gravity, Leukocytes, Ascorbic Acid



WATERBATH 10L

- One-hole
- Power: 300w
- Volume :10L