



Orthopaedic



ABS ORTHOPAEDIC TRACTION BED

- Dimension: L2020xW900xH500mm
- Bed base is made of cold-rolled steel plate and coated by electrostatic spray.
- Aluminum alloy guard rails, retractable, easy operation.
- Built-in cranks control back rest & leg rest upward and height adjustable function.
- With protection rubber at the leg end. With anti-bump cushion device on the bed corners for protection against to wall or etc..
- With multi-function orthopedic traction frame, ensure traction treatment.



AMBULANCE STRETCHER

- Product Size (L x W x H): High position: 190 x 55 x 92cm; Low position: 190 x 55 x 26cm
- Maximum Angle of the back: 70°
- Carton size (1pc/carton): 195 x 62 x 28cm
- N.W.: 36.8kg GW.: 43kg
- Load bearing: less than 159kg"



AXILLIARY/ARMPIT CLUTCHES SIZE S

- Mobility Aid for Non-weight-bearing recovery
- Lightweight Aluminum construction.
- Adjustable underarm pad and hand grip.
- Replaceable PU foam pad.
- Size: small

To Request for Quote:

WhatsApp: +254 739 643897



CEREBRAL PALSY WHEELCHAIR

- Wheelchair width: 56cm
- Folding width: 33cm
- Seat width: 33cm, depth: 39cm, height 48cm
- Wheelchair height 98cm, length 106cm
- Rear wheel diameter 16", front wheel diameter 6"
- Backrest height 40cm
- Max capacity 100kg
- Weight 22.86kg



CLIMBING CHAIR

- With track, only one person can push down the stretcher telescoping backrest.
- It is made of high strength powder-coated aluminum alloy PVC seat.
- Unfold size: 80*66*163cm
- Folded size: 106*66*34cm
- Packing size: 110*70*38cm
- G.W: 45kg
- Load bearing: 160kg



HINGED KNEE BRACES OPEN PATELLA

- Material: fabric neoprene.
- Size: universal
- Usage: for knee support



WALKING FRAME WITH REMOVABLE WHEELS

- Lightweight Folding Wheeled Walking Frame
- Lightweight Aluminium Construction
- Adjustable Height
- Easy-to-Clean Plastic Handgrips
- Non-Marking, Slip-Resistant Rubber Tips