



Radiology Equipment



AUTOMATIC X-RAY FILM PROCESSOR

- 2.4 inch LCD screen, can display developer, fixer, drying temperature and processor speed.
- With Intel chip, it has the function of failure analysis and troubleshooting.



ECG MACHINE

- Sync collection for 12-lead ECG, adopt digital 1 signal processing technology and get high- quality ECG waveform via power frequency filter, baseline filter and EMG filter of ECG signal.
- Display of 3/6/12-lead ECG, print mode, sensitivity, paper speed and filter state, etc. on one screen, convenient for contrastively diagnosing."



LAPTOP ULTRASOUND SCANNER

- Pretty and lightweight
- 15 inch LCD screen
- Full digital imaging technology, crystal-clear images
- Broadband multi-frequency probes
- 2 probe connector
- Silica gel backlight gravure laser keyboard
- ARM system powerful function, include report, can work with laser printer.
- Built-in battery (can work more than 3 hours)

To Request for Quote:

WhatsApp: +254 739 643897



LEAD PROTECTIVE GLOVES FOR X-RAY

- Radiation protection: Offers protection against harmful X-rays or other ionizing radiation.
- Flexibility and dexterity: Designed to allow for movement and precision while handling instruments during medical procedures.
- Comfort: Padded interiors to ensure comfort during extended use.



MOBILE X-RAY MACHINE

- Adopting high frequency inverting technology, with stable high voltage output to ensure quality images.
- Compact design makes it easy and portable for usage in different locations.
- With 3 modes of exposure control: Remote, hand switch and panel key.
- Break self-diagnostic and self-protecting
- Can be upgraded to digital x-ray system (YSX050-C) simply by adding a flat panel detector.



PORTABLE ULTRASOUND

- Strong and durable ultrasound in a portable carry case, gel and an AC adapter.
- Preheating function for comfort
- Pulsed and continuous modes
- 3 output intensities
- 3 timer settings



PRINT PAPER 110MM*18M

- Paper size: 110 mm (W) x 18 m (A6 width).
- Print quantity: 250 prints/roll.
- Number of rolls or pack: 10 rolls.
- Water resistance and high-gloss: The high gloss/gloss layer of print media, provides high water resistance and high storage stability.
- Superior print quality: Thanks to rigorous application pressure control, the thermal coat layer delivers high-quality coloring properties. The curve and Dmax are strictly adjusted to ensure the stable provision of consistent, optimal image quality.



SONY ULTRASOUND PRINTER

- Input Current: 0.6Amp
- Frequency: 50 Hz
- Dimensions (W x H x D): 154 x 88 x 240 mm (6 1/6 x 3 1/2 x 9 1/2 inches)
- Mass: 2.5 kg (5 lb 8 oz)
- Model No: UP-X898MD
- Brand: Sony
- Voltage: 240 V
- Printer Driver Software: Microsoft Windows 7 32/64 bit, 8 32/64 bit, 10 32/64 bit Gradations: 8 bits (256 levels) processing
- Resolution: 325 dp
- Picture Elements: Digital: 4096 x 1280 dots, Video NTSC: 720 x 504 dots, Video PAL: 720 x 604 dots



X-RAY PROTECTIVE CLOTH

- "0.35mmPb
- Size: 1100*620mm, 1000*600mm, 900*580mm"
- Double side with half sleeve



XRAY VIEWER-SINGLE

- "For 2 units film (85*42cm) 120W
- LED light source
- Electrophoretic aluminum alloy frame
- High definition light penetration
- Acrylic Screen
- Easy maintenance
- Inside clipper device, film clipped firmly and taking out easily."



X RAY PROTECTIVE COLLAR LEAD THYROID COLLAR

Features and Benefits:

- Thyroid protection during , dental x-ray and mammography procedures.
- Magic Stick adjustable, one size fits all.
- Light-weight and comfortable.
- Personal thyroid protection.
- 0.50mm PB meets current FDA ASTM protective standards