



Theatre Equipment



ANESTHESIA MACHINE

- 10.4 inch screen
- Ventilation Mode: IPPV, SIGH, SIMV/F/2, SIPPV (Optional: PEEP), MANU
- O₂/N₂O flow tube flow meter
- Tidal volume: 30-150ml.
- Two Gas Supplies: O₂/N₂O: 0.3MPa-0.5MPa
- Vaporizer Quantity: Isoflurane, Sevoflurane, Halothane (Optional two of three)
- Application: Adult, children



BIPHASIC DEFIBRILLATOR

- Non-synchronizer: Outlife defibrillator
- Sine wave: Biphasic technology
- Energy: 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 joules
- Charging time: at 360J less than 10 sec
- Paddle options: reusable external adult paddles (pediatric paddle integrated)
- Battery: operates from the mains (AC) and with its internal Lead-acid batter (Li-ion optional)



DIATHERMY MACHINE 400W

- The output power: 400W
- Six output modes: Cutting (400W) Blend cut1 (200W), Blend cut2 (150W), Soft Coag (120W), Strong Coag (100W), Bipolar Coag (70w), Atmospheric pressure range 86.0-106.0 kPa
- Power supply: 220V±22V, 50Hz±
- Working frequency: 360-460kHz
- The rated power of the equipment: 880VA+10%
- Type of device: CF
- The machine power consumption: ≤100VA

To Request for Quote:

WhatsApp: +254 739 643897



DRESSING TROLLEY

- Metallic body casing,
- It has four 75mm castors
- Has a drawer
- Length 23 inch
- Width 15 inch
- Height 19.5 inch
- Bowl 73 inch



ECG 12 CHANNELS

- Electrocardiograph 12 channel
- It is a machine that shows and record heart rhythm and electrical activity.



ELECTRIC THEATRE OPERATING TABLE

- Table top lifting
- Forward and backward tilting
- Leftward and rightward tilting
- With emergency power



HAND HELD PULSE OXIMETER DISPLAYS SP02 AND PULSE WAVEFORMS

- Plethysmograph graph displayed for improved performance
- Color TFT Display with Multi-Directional Views
- Brightness control
- Low power consumption with low battery alarm on display.
- Adult SpO2 sensor Included.
- Optional sensors for neonatal and pediatric applications
- Data analysis software
- Programmable alarms
- AC charger or from USB port
- Programmable display
- Flash memory with 24 data read back



LARYNGOSCOPE

- Size: 2#, 3#, 4#
- Basic assembly (per 1 set)
- No of blades: 3 Pcs
- Handles: 1 PC
- Bulb: 3PCS (2.5V)



MAYO TROLLEY

- Medical tray: tray with a smooth surface and durable.
- Durable frame: thickened stainless steel pipe welding, more durable.
- Adjustable knob: height adjustable by rotation knob.
- Stable base: thickened ss square tube more stability.
- Silent universal castor: pressed steel housing with pu cover water resistance.



OPERATING LIGHT HALOGEN STAND/ CEILING

- Theatre light also known as surgical light head or theatre light.
- Medical device intended to assist medical personnel during a surgical procedure.



OXYGEN FLOMETER (VERTICAL CONNECTION)

- Material: copper mainbody, plastic handwheel
- Input Pressure: 150bar
- Output Pressure: 2 bar-3bar
- Control pressure of safety valve: 3.5bar $\pm$ 0.5bar
- Flowmeter control range: 0-10LPM/15LPM
- Material: Main brass + Aluminum alloy flow meter
- Connecting thread: Male G5/8 (top-thread)